



Electrotherapy & Physical Modalities

PRODUCT CATALOGUE



Scan to Download
the Catalogue

YOUR RELIABLE & TRUSTED PARTNER

**In Physiotherapy, Rehabilitation
& Occupational Health Solutions**



ABOUT US

IPEX Healthcare Limited is a Tanzanian based and Tanzanian registered company dealing with Importation, and Supply of Hospital consumables, Medical, Surgical and Laboratory equipment, supplies, devices and other related products and services in the medical and healthcare sector.

LGT-2500X

Electrical Extracorporeal
Shock Wave Therapy (ESWT)
Device LONGEST



★ FEATURES

- Modern and sleek appearance design
- Sturdy metal shell, durable structure
- Integrated with vibration therapy, frequency up to 35Hz
- Control panel integrated into the shockwave handpiece
- A wide range of preset protocols with image guidance
- User-defined protocols for personalized treatment
- Automatic recognition of connected handpiece
- Patient management system with treatment history
- Record the number of shocks delivered by handpieces

★ HIGHLIGHTS

- Large touchscreen
- Rotary knob for quick parameters setting
- Energy range from 0.3 bar
- Automatic handpiece recognition
- Handpiece with a control panel
- Preset protocols with image guidance
- User-defined protocols
- VAS & ROM assessment system
- Patient database
- Account management to safeguard patient data



LGT-2520GP


Electromagnetic Focused Shockwave Therapy Device LONGEST



LGT - 2520GP utilizes the advanced "Electromagnetic Shockwave" technology, which generates a pressure field that converges in the adjustable focus at determined depths in specific tissues where maximal pressure is achieved. It can accurately target the root cause of musculoskeletal conditions, effectively relieve pain, promote tissue healing, and reduce inflammation.

★ HIGHLIGHTS

- Electromagnetic shockwaves for accurate treatment
- Easy adjustment of shockwave energy level and frequency
- Soft, skin-friendly coupling area, closely contact focal area
- Adjustable coupling area to precisely treat the targeted spot at different depths
- An inbuilt protocol library equipped with a 3D body image helps locate the treatment site
- A patient database creates and stores individualized treatment plans and treatment histories for each patient



★ FEATURES

- Electromagnetic shockwaves for accurate treatment
- Frequency up to 20Hz
- Easy adjustment of shockwave energy level and frequency
- Soft, skin-friendly coupling area, closely contact focal area
- Adjustable coupling area to precisely treat the targeted spot at different depths

LGP-2500 SERIES

Electrical Extracorporeal Shock Wave Therapy Device LONGEST POWER SHOCKER – LGP-2500S & LGP-2500S Plus longest



LGT-2500S

- 7" touchscreen with an intuitive interface
- User-defined protocols
- Pressure: 1.0-4.0 bar
- Frequency: 1-17 Hz
- Pulse: 1-3,000



LGT-2500S Plus

- Integrated vibration therapy
- Pressure: 1.0-5.0 bar
- Frequency: 1-22 Hz
- Pulse: 1-8,000

★ APPLICATIONS

- Anti-cellulite treatment
- Body reshaping
- Skin elasticity improvement
- Treatment of skin irregularity after liposuction
- Connective tissue tightening
- Scar and wrinkle smoothing
- Treatment of stretch marks
- Lymphatic system stimulation



★ TECHNICAL SPECIFICATIONS

Model Number	LGT-2500S PLUS	LGT-2500S
Image		
Channel	2 1 for Shockwave handpiece, 1 for Vibration handpiece	1
Display Screen	7" LCD Touch	7" LCD Touch
Pressure Range	1.0~5.0 Bar	1.0 4.0 Bar
Frequency	1~22Hz	1~17Hz
Vibration Handpiece:	1~35Hz	/
Pre-set Protocols	Yes	Yes
User-defined Protocol	Yes	Yes
Item Dimension (Without handpiece holders)	39.5 (L) x 39 (W) x 19.4 (H) cm	39.5 (L) x 39 (W) x 19.4 (H) cm
Item Weight (Handpiece Excluded)	11.3kg	10.5kg
Languages	10	7

LGT-3500 Series

High Frequency therapy Device LONGEST

LGT-3500 high frequency therapy device redefines non-invasive electrotherapy. Powered by advanced TECAR Therapy technology, it delivers targeted, safe, and effective treatment for pain management, tissue repair and more—tailored to meet the diverse needs of clinics and hospitals.



LGT-3500LP



LGT-3500GP



★ LGT-3500GP & LGT-3500LP HIGHLIGHTS

- **Protocol Management**
Picture guidance (treatment area + operation method), and customized prescription function.
- **Customizable Output**
Adjustable intensity & duration for personalization.
- **Real-Time Monitoring**
Temperature, resistance value, voltage value, power.
- **Smart Safety**
Auto-identifies accessories (reduces burn risk) + over-temperature alarm+open-circuit alarm
- **Patient Record Management**
Stores history, supports multi-account, and protects privacy.
- **Multiple Accessories**
High-quality accessories are optional to adapt to diverse clinical needs.





LGT-2210DS and LGT-2200SP

Compression Therapy Devices LONGEST



★ HIGHLIGHTS

- Smart garment recognition
- Turn off chambers individually
- Large touchscreen with an intuitive user interface 29 preset compression modes
- Wide range of overlapping cuffs



LGT-2200SP Compression Therapy Device

★ HIGHLIGHTS

- Portable & battery operated
- 6 compression modes
- Full-body compression experience

LGT-2320SP

Electrotherapy Stimulation Therapy Device (Portable) LONGEST



★ HIGHLIGHTS

- 7-inch touch screen
- 7 preset protocols with electrode placement guidance and treatment guidance
- Custom protocols for personalized treatment
- "12+1" channels for improved efficiency
- Convenient self-adhesive electrode pads of various sizes Various types of pen electrodes
- Patient database with treatment records
- Quiet operation



★ APPLICATIONS

▪ Muscle Warm-up Aid

The Electro Stimulation Therapy Device proves to be highly beneficial for pre-workout warm-ups. The Electro Stimulation Therapy Device sends gentle electrical impulses to muscles, enhancing blood circulation and oxygen supply to muscle tissues. This loosens and readies the muscles for intense activity and activates the neuromuscular system to improve nerve-muscle communication, muscle reactivity, and coordination.

▪ Muscle Training & Strengthening

The Electro Stimulation Therapy Device can conduct targeted training for specific muscle groups. The intensity and frequency of electrical impulses can be adjusted to suit different fitness levels and training goals, and by strengthening muscles and improving endurance, it can enhance athletic performance in various sports and activities.

▪ Post-exercise Muscle Relief

It offers significant relief for post-exercise muscle soreness. The device sends gentle electrical impulses to fatigued muscles, increasing blood flow to speed up lactic acid and waste removal. It boosts the body's recovery, reducing inflammation and muscle stiffness. Athletes and fitness lovers can use it for quicker post-exercise recovery, shortening downtime, and maintaining fitness and performance.

▪ Injury Rehabilitation

The device helps in the recovery process by stimulating the affected muscles and promoting blood circulation. It can reduce pain and inflammation, accelerate healing, and restore muscle strength and flexibility. Additionally, it aids in breaking up muscle knots and adhesions, further facilitating the return to optimal muscle function and range of motion.





LGT-2320ME

**Electrotherapy Stimulation
Therapy Device (Portable)
LONGEST**

A professional low-frequency electrical stimulation station for all stages of neuromuscular dysfunction and dysphagia.

★ HIGHLIGHTS

- 7-inch touch screen
- 31 preset protocols with electrode placement guidance and treatment guidance
- Custom protocols for personalized treatment
- "12+1" channels for improved efficiency
- Convenient self-adhesive electrode pads of various sizes
- Various types of pen electrodes
- Patient database with treatment records
- Quiet operation



LGT-233

Portable Electro-Stimulation Therapy Device LONGEST



★ HIGHLIGHTS

- Lightweight & Wearable
- User App Control "MStim Drop"
- Built-in Sensor to Monitor User Movement and Provide Real-time Feedback Certification: FDA, Health Canada

★ FEATURES

▪ Built- in Smart Sensor

The built- in gyroscope and acceleration sensor control the timing and duration of electrical stimulation by tracking the swing angle and pace of patient's leg.

▪ Wireless Bluetooth Connection

Quick connection easy and convenient

▪ Small Unit

The unit is small and portable, and use with straps. It can be used for a long time without any sense of heavy while walking.

▪ Real- time Dynamic Observation of Walking Angle

Observe walking angle changes every moment and record walking distance and number of steps simultaneously to facilitate quantitative motion data.

▪ User- defined Programs

Customized protocols for per patients based on their conditions

TECHNICAL SPECIFICATIONS

Mode	Train Mode, Walk Mode
Output Waveform	Symmetrical biphasic pulse
Pulse Duration	50- 500µs
Pulse Frequency	1- 120Hz
Intensity	0- 100mA
Battery	Rechargeable Lithium Battery




LGT-231

Portable Muscle Stimulator TENS & NMES Device LONGEST



Neuropathic Pain



Dysphagia

★ HIGHLIGHTS

- Compact and Portable
- Battery-Operated
- 2 Channels
- NMES & TENS
- Neuropathic Pain
- A Wide Range of Indications
- Preset Protocols with Pad Placement Guidance Custom Protocols
- User App Control - "Mstim Reha"
- Certification: CE, FDA, Health Canada
- 2 Years of Warranty

★ FEATURES


- **Wireless Bluetooth connection**
Quick connection, easy and convenient
- **User-defined Programs**
Parameters such as frequency, pulse width, treatment time are adjustable in all treatment modes
- **Preset protocols**
Built in parameters and schematic diagram of muscle patch

TECHNICAL SPECIFICATIONS

	CE Version	FDA Version 510(K) Number: K182020
Mode	P1/P2/P3	TENS/NMES
Pulse Duration	50- 400Qs	50- 400Qs
Pulse Frequency	2- 150Hz	1- 120Hz
Intensity	1/P2: 0- 60mA P3:0- 1000QA	0- 100mA
Time	1- 60min	1- 60min
Battery Life	More than 4 hours	More than 4 hours

MSTIM REHA LGT-231

Portable Electro-Stimulation Therapy Device Mstim Reha LGT-231 is a lightweight and portable multifunctional electrotherapy device that provides two channels of TENS, NMES current to stimulate nerves or muscles to help relieve pain and improve neuromuscular function





LGT-2900EUV

Combination Therapy
Electrotherapy, Ultrasound,
& Vacuum Therapy LONGEST

LGT-2900EUV

Ultrasound and Electrotherapy Device

★ HIGHLIGHTS

- Intuitive 10.1" full-color touchscreen interface
- Modern and sleek appearance design
- Sturdy metal shell, durable structure
- 3 treatment modalities combination
- Portable, lightweight, battery-operated ultrasound module, used individually
- Rotary knob for quick parameters setting
- Easy to move with an optional trolley
- A wide range of preset protocols
- User-defined protocols for personalized treatment Patient management system with treatment history



LGT-2900U

Ultrasound Therapy Device



★ FEATURES

- Compact dimensions, 2.8kg battery-operated, good portability
- Built-in rechargeable battery (9000mAh)
- Intuitive 10.1" touchscreen interface
- 42 preset ultrasound protocols with body part navigation
- User-defined protocols
- 1cm² and 5 cm² applicators, lightweight & ergonomic
- Patient database and account management system
- Multi-language software
- Software update via USB
- *Carrying bag
- *Trolley

LGT-2900EU

Ultrasound Electrotherapy Device



★ FEATURES

- 2 treatment modalities: ultrasound and electrotherapy
- Ultrasound and electrotherapy combination treatment
- 4-channel electrotherapy, treat multiple people or body parts simultaneously
- Compact and battery-operated, great portability
- Built-in rechargeable battery (9000mAh)
- Intuitive 10.1" touchscreen interface
- 42 preset protocols for ultrasound, 57 for electrotherapy
- User-defined protocols
- 1cm² and 5 cm² applicators, comfortable and easy to handle
- A wide range of advanced current waveforms available
- Self-adhesive electrodes & silicon rubber electrodes of various sizes
- Patient database and account management system
- Software update via USB
- *Carrying bag
- *Trolley

LGT-2900EUV

Ultrasound Electrotherapy Device with Vacuum Unit



★ FEATURES

- 3 treatment modalities: ultrasound, electrotherapy, and vacuum therapy
- Multiple channels, treat multiple people or body parts simultaneously
- Ultrasound and electrotherapy/vacuum therapy combination treatment
- Compact and battery-operated, great portability
- Built-in rechargeable battery (9000mAh)
- Intuitive 10.1" touchscreen interface
- 42 preset protocols for ultrasound, 57 for electrotherapy, 25 for combination
- User-defined protocols
- 1cm² and 5 cm² applicators, comfortable and easy to handle
- A wide range of advanced current waveforms available
- Self-adhesive electrodes & silicon rubber electrodes of various sizes
- High-quality suction cups (vacuum electrodes) of various sizes
- Patient database and account management system
- Software update via USB
- *Carrying bag
- *Trolley



TECHNICAL SPECIFICATIONS

Model		LGT-2900U	LGT-2900EU	LGT-2900EUV
Ultrasound	Channel	2-channel		
	Acoustic Working Frequency	1 MHz and 3 MHz		
	Waveform	Continuous wave, pulsed wave		
	Pulse Frequency	16 Hz, 48 Hz and 100 Hz		
	Intensity	In Continuous mode: 0-2 W/cm² (1 cm² Applicator) 0-2.5 W/cm² (5 cm² Applicator) In Pulsed mode: 0-3 W/cm² (1 cm² Applicator) 0-5 W/cm² (5 cm² Applicator)		
	Treatment Time	1-30 minutes		
Electrotherapy	Channel	—	4-channel	
	Output Mode	—	Constant Current (CC) or Constant Voltage (CV)	
	Treatment Time	—	1-60 minutes	
Vacuum Therapy	Channel	—	—	4-channel
	Pressure Range	—	—	10-40 kPa
	Vacuum Mode	—	—	15 times per min/continuous
Dimensions		379×310×146 mm	379×310×211 mm	379×310×261 mm
Weight		2.8kg	4.1kg	5.6 kg
Interface		10.1" touch screen		
Power Supply		Battery: 11.1V, 9000mAh, lithium AC 100-240V, 50/60 Hz		

LGT-2410S


Cryotherapy Device LONGEST

★ HIGHLIGHTS

- 10.1-inch full-color LCD touchscreen
- Therapeutic encyclopedia with image guidance
- Real-time detection of skin temperature
- Hypothermia pop-up with automatic stop and light cues
- NRS function to track the treatment progress
- Integrated trolley for better mobility

TECHNICAL SPECIFICATIONS

Interface	10.1" color touchscreen
Treatment Duration	1s-10min
Dimensions	40(L) x 48(W) x 125(H)cm
Weight (without cylinder)	29kg
Gas Cylinder Capacity	6/8/10L
Applicator	● Screw thread G5/8 ● Monitoring screen ● Temperature sensor ● LED light ● Anti-frosting nozzle
Gas Cylinder Requirement	● External Diameter < 175mm ● Height < 870mm ● with a Rising Tube/Siphon ● Food-grade CO2
Functions	● Real-time skin temperature monitoring ● Treatment image guidance ● Operation guidance ● Automatic shut off ● Patient database ● Data export via USB ● Treated area positioning ● NRS function ● Cylinder capacity monitoring
Temperature	0 °C-15 °C
Temperature Detection	-3 °C-40 °C
Power Supply	AC100-240V, 50/60Hz




LGT-2360S

Electrostatic Oscillation Therapy LONGEST PowerOsci

★ HIGHLIGHTS

- Handheld and rechargeable
- Preset protocols
- Easy to learn, easy to use
- Low maintenance
- Suitable for both healthcare professionals and individuals
- FDA-cleared



TECHNICAL SPECIFICATION

Discover how you can stand out from the crowd in electrostatic oscillation therapy with PowerOsci LGT-2360S

Carrying Bag

convenient to take when going out

2.8" Multiple Color Display

for stunning visuals and better readability

10cm Applicator

to treat larger areas easily, improving efficiency

5.5cm Applicator

with easy-to-replace head for hygiene and patient safety

Applicator Handle

ergonomic design, comfortable to hold, agile to operate

Device Stand

to top the device on for easy charging and storage

Electrode Pads

adhesive, easy to remove, no residue

Size: 76.5 x 143.5 x 35.5mm

Weight: 230g

Output Frequency: 5-300 Hz

Output Voltage: Max: 450V

Rechargeable Battery: 2600mAh, Last for 6h, Recharged in Less than 3h

Quality Assurance: FDA approved; 2-year warranty





LGT-2202DVT

VTE Prevention Pump LONGEST

LGT-2202DVT is a state-of-the-art medical device designed to provide highly effective and targeted compression therapy for DVT patients. The device is remarkably compact in dimensions and easy to operate, highly suitable for hospitals to simplify workflow and enhance patient outcomes.



★ HIGHLIGHTS

- Portable & Lightweight Battery-operated
- Quiet operation
- Automatic garment recognition
- Alternate and synchronous inflation
- One leg compression supported
- Vascular refill detection, safe and smart compression cycles
- Real-time pressure monitoring
- Various kinds of sleeves (single-patient use)

SPECIFICATION

Chambers	1-3
Mode	4
Time Range	Continuous
Pressure Range	30 - 130mmHg
Interface	2.8" LCD screen
Rechargeable Battery	Yes (Lithium battery: 3,500mAh)
Weight	1.5kg
Dimension	145(L) x 166(W) x 188(H)mm
Single-leg Treatment	Supportive
Sleeve Auto-Detection	Supportive
Vascular Refill Detection	Supportive

★ FEATURES

- Compact dimensions and lightweight
- Battery-operated
- Quiet operation
- Automated garment recognition Intuitive user interface
- Alternate and synchronous inflation Single leg treatment supported
- Blood velocity detection, safe and smart compression cycles
- Real-time pressure monitoring
- Gradient compression
- High-quality connector, easy to connect and disconnect
- Complete range of garments, made of high quality and breathable fabrics



One-Button-Start System

Automatically setting the proper pressure and compression cycle, reducing additional user intervention.



Smart Sleeves Recognition

Auto-Recognize the using sleeves belonging to Longest DVT sleeves family.



Compatible Bedrail Hanger

Lightweight, battery-operated. Increase the general mobility and enabling medical staff to move it between patients from different bed types.



Quick and Safe Connection

Equipped with a convenient connector to ensure the tightness of the sleeve connection.



Noiseless Operation

Ensure a quiet environment (below 42dB) for patients' rest.



Vascular Refill Detection

Detect each patient's vascular refill status to deliver a customized compression therapy plan.

BENEFITS OF THE TREATMENT



VTE
Prevention



Circulation
Enhancement



Acute and Chronic
Edema Reduction



Venous Stasis
Treatment



Improve the Healing
Speed of Leg Ulcers



Venous Insufficiency
Prevention and Treatment



LGT-2201DVT & LGT-2200DVT

VTE Prevention Pump –
Compression Therapy
Device LONGEST

★ FEATURES

- 4 Compression Modes
- 5" Touch Screen
- Reusable & Disposable Sleeves
- Gradient Pressure Setting
- Intermittent Pneumatic Compression >
- Pressure of Every Chamber Adjustable
- Rechargeable Battery Powered by AC or DC

TECHNICAL SPECIFICATIONS

Chambers	3
Compression Mode	4
Time Range	1min-24hours or continuous
Pressure Range	40-150mmHg
Hold Time	0-10s(Step by 1s)
Interval Time	10-60s(Step by 1s)
Rechargeable Battery	Yes (Lithium battery : 11.1V 9000mAh) Maximum Running Time: No Less than 8 Hours
Adaper	Yes
Input Voltage	AC100-240V 50-60Hz
Noise Level	<50 dB
Dimension	250(L)x184(W)x150(H)mm
Weight	2.5kg

LGT-2200DVT

Compression Therapy Device

- Rechargeable battery (9000mAh)
- 4 compression modes
- Wide pressure range: 40-150 mmHg
- Screen lock to prevent mistaken touching



LGT-2201DVT

Compression Therapy Device

- Easy button controls
- Rechargeable battery (3500mAh)
- 4 compression modes
- Wide Pressure Range: 30-150 mmHg



VTE

VTE preventive sleeve LONGEST




Order No.	Sleeve	Size	Proximal Circumference
220092	Thigh Length Sleeve	S	<54.9 cm
220092M		M	<69.5 cm
220092L		L	<97.3 cm
220093M	Calf and Foot Sleeve	M	<57.0 cm
220093L		L	<72.9 cm
220094	Calf Sleeve	S	<61.0 cm
220094M		M	<71.0 cm
220094L		L	<80.1 cm
220095	Foot Sleeve	S	< 45.2cm
220095M		M	< 55.6cm
220095L		L	<69.3cm
220098M		M	<38.4 cm

Order No.	Sleeve	Size	Proximal Circumference
220292	Thigh Length Sleeve	S	<54.9 cm
220292M		M	<69.5 cm
220292L		L	<97.3 cm
220293	Calf and Foot Sleeve	S	<54.0 cm
220293M		M	<57.0 cm
220293L		L	<72.9 cm
220294	Calf Sleeve	S	<61.0 cm
220294M		M	<71.0 cm
220294L		L	<80.1 cm
220298M	Foot Sleeve	M	<38.4 cm

LGT-2660 Series

Electromagnetic Stimulator – Musculoskeletal Pain

LGT-2660 Series Electromagnetic Stimulator for Musculoskeletal Pain is an advanced rehabilitation device engineered to address a wide range of musculoskeletal pain conditions and associated dysfunctions. Operating on the principle of electromagnetic induction, It generates a magnetic field and inductive currents to regulate neuromuscular function while delivering deeply penetrating energy—this energy transcends superficial tissues to act directly on muscles, joints, and bones. For both acute and chronic musculoskeletal pain, it provides a safe, efficient, and medication-free new solution.



★ APPLICATIONS AREAS

Musculoskeletal & Neurological Rehabilitation

- Hospitals (Orthopedics, Rehab, Pain Management)
- Clinics

Sports Rehabilitation

- PT Centers
- Rehab Centers
- Sports Medicine Facilities

Geriatric Rehabilitation

- Hospitals
- Clinics
- Rehab Centers
- Nursing Homes

LGT-2660AP

Electromagnetic Stimulator (Single-channel)



★ OPTIONAL MODELS

- 15.6-inch full-color touch screen with a knob for quick parameter adjustment
- Equipped with a universal stand for multi-angle adjustment, adapting to diverse patient requirements
- Picture guidance to assist in locating stimulation areas
- Adjustable treatment parameters to meet different therapeutic needs
- Preset and custom prescription functions for flexible therapy adaptation
- Supports binding patient information and saving treatment data
- Equipped with a trolley for easy mobility

LGT-2660BP

Electromagnetic Stimulator (Dual-channel)





LGT-2640AP

Electromagnetic Stimulator LONGEST – Pelvic Floor and Urinary System Disorder

LGT-2640AP Magnetic Stimulator specializes in treating pelvic floor and urinary system disorders. Its non-invasive magnetic fields target pelvic muscles and nerves, accelerating symptom relief and recovery.

FEATURES

Intuitive 15.6" color touchscreen


Ergonomic seat featuring a built-in coil and easy installation

24 preset protocols with anatomical diagrams

Auto-save for case data & treatment parameters

Instant settings with Rotary Knob

Emergency stop button for safety




Indications for Women

Primary Causes

- Pregnancy, Childbirth, Aging

Key Conditions

- Stress Urinary Incontinence (SUI)
- Urgency Urinary Incontinence Overactive Bladder (OAB)
- Pelvic Organ Prolapse (POP)



Indications for Men

Primary Causes

- Prostate Disease, Surgical Injury, Aging

Key Conditions

- Stress Urinary Incontinence (Postoperative)
- Chronic Prostatitis
- Chronic Pelvic Pain Syndrome (CP/CPPS)
- Erectile Dysfunction (ED) – Adjuvant Therapy


LGT-2410S

Cryotherapy Device LONGEST

Freeze the Pain Away! The cryotherapy device LGT-2410S uses advanced localized cryotherapy technology, offering unmatched precision, efficacy, safety, and ease of use. The device combines the effects of cold and pressure.

★ HIGHLIGHTS

- 10.1-inch full-color LCD touchscreen
- Therapeutic encyclopedia with image guidance
- Real-time detection of skin temperature
- Hypothermia pop-up with automatic stop and light cues
- NRS function to track the treatment progress
- Integrated trolley for better mobility




TECHNICAL SPECIFICATIONS

Interface	10.1" color touchscreen
Treatment Duration	1s-10min
Dimensions	40(L) x 48(W) x 125(H)cm
Weight (without cylinder)	29kg
Gas Cylinder Capacity	6/8/10L
Applicator	▪ Screw thread G5/8 ▪ Monitoring screen ▪ Temperature sensor ▪ LED light ▪ Anti-frosting nozzle
Gas Cylinder Requirement	▪ External Diameter<175mm ▪ Height<870mm ▪ with a Rising Tube/Siphon ▪ Food-grade CO2
Functions	▪ Real-time skin temperature monitoring ▪ Treatment image guidance ▪ Operation guidance ▪ Automatic shut off ▪ Patient database ▪ Data export via USB ▪ Treated area positioning □NRS function ▪ Cylinder capacity monitoring
Temperature	0 °C-15 °C
Temperature Detection	-3 °C-40 °C
Power Supply	AC100-240V, 50/60Hz

LGT-2320SP

Electrotherapy Stimulation Therapy Device (Portable) LONGEST



★ HIGHLIGHTS

- 7-inch touch screen
- 7 preset protocols with electrode placement guidance and treatment guidance
- Custom protocols for personalized treatment
- "12+1" channels for improved efficiency
- Convenient self-adhesive electrode pads of various sizes
- Various types of pen electrodes
- Patient database with treatment records
- Quiet operation



7" Full-color Touch Screen

Lightweight and Portable

New UI Design with Easy User Interaction

Rotary Knob for Quick Adjustment

Image Guidance of Electrode Placement

13 Channels for Multisite Treatment

Quiet During Operation

Pen Electrode as a Standard Configuration

An Emergency Switch for Treatment Safety

LGT-2320ME

Electrotherapy Stimulation Therapy Device (Portable) LONGEST

A professional low-frequency electrical stimulation station for all stages of neuromuscular dysfunction and dysphagia.



★ HIGHLIGHTS

- 7-inch touch screen
- 31 preset protocols with electrode placement guidance and treatment guidance Custom protocols for personalized treatment
- "12+1" channels for improved efficiency
- Convenient self-adhesive electrode pads of various sizes
- Various types of pen electrodes
- Patient database with treatment records
- Quiet operation



FEATURES

Intuitive 15.6" color touchscreen

Ergonomic seat featuring a built-in coil and easy installation

24 preset protocols with anatomical diagrams

Auto-save for case data & treatment parameters

Instant settings with Rotary Knob

Emergency stop button for safety



LGT-2520GP

Electromagnetic Focused Shockwave Therapy Device LONGEST



LGT-2520GP utilizes the advanced “Electromagnetic Shockwave” technology, which generates a pressure field that converges in the adjustable focus at determined depths in specific tissues where maximal pressure is achieved. It can accurately target the root cause of musculoskeletal conditions, effectively relieve pain, promote tissue healing, and reduce inflammation.

★ HIGHLIGHTS

- Electromagnetic shockwaves for accurate treatment
- Easy adjustment of shockwave energy level and frequency
- Soft, skin-friendly coupling area, closely contact focal area
- Adjustable coupling area to precisely treat the targeted spot at different depths
- An inbuilt protocol library equipped with a 3D body image helps locate the treatment site
- A patient database creates and stores individualized treatment plans and treatment histories for each patient

★ FEATURES

- Electromagnetic shockwaves for accurate treatment
- Frequency up to 20Hz
- Easy adjustment of shockwave energy level and frequency
- Soft, skin-friendly coupling area, closely contact focal area
- Adjustable coupling area to precisely treat the targeted spot at different depths

LGP-2500 SERIES

Electrical Extracorporeal Shock Wave Therapy Device LONGEST POWER SHOCKER - LGP-2500S & LGP-2500S Plus



PowerShocker LGT-2500S

Electrical Extracorporeal ShockWave Therapy Device

★ HIGHLIGHTS

- 7" touchscreen with an intuitive interface
- User-defined protocols
- Pressure: 1.0-4.0 bar
- Frequency: 1-17 Hz
- Pulse: 1-3,000



PowerShocker LGT-2500S Plus

Electrical Extracorporeal
ShockWave Therapy Device

★ HIGHLIGHTS

- Integrated vibration therapy
- Pressure: 1.0–5.0 bar
- Frequency: 1–22 Hz
- Pulse: 1–8,000

TECHNICAL SPECIFICATIONS

Model Number	LGT-2500S PLUS	LGT-2500S
Image		
Channel	2 1 for Shockwave handpiece, 1 for Vibration handpiece	1
Display Screen	7" LCD Touch	7" LCD Touch
Pressure Range	1.0–5.0 Bar	1.0 – 4.0 Bar
Frequency	1–22Hz	1 – 17Hz
Vibration Handpiece:	1–35Hz	/
Pre-set Protocols	Yes	Yes
User-defined Protocol	Yes	Yes
Item Dimension (without handpiece holders)	39.5 (L) x 39 (W) x 19.4 (H) cm	39.5 (L) x 39 (W) x 19.4 (H) cm
Item Weight (Handpiece Excluded)	11.3 kg	10.5 kg
Langages	10	7



PowerWave LGT-2500S

The ideal body shaping treatment
Acoustic waves provide an ideal body shaping
treatment option

TECHNICAL SPECIFICATIONS

Pressure Range: 1.0–4.0 Bar (Step by 0.1 Bar)

Frequency: 1–17Hz

Pre-set Protocol: Yes

Circumference: Yes

Dimension: 39.5 (L) x 39 (W) x 19.4 (H) cm

Standard Weight: 11.3kg

Langages: English, Italian, German, Spanish, French

LGT-3500 Series

High Frequency therapy Device

LGT-3500 high frequency therapy device redefines non-invasive electrotherapy. Powered by advanced TECAR Therapy technology, it delivers targeted, safe, and effective treatment for pain management, tissue repair and more—tailored to meet the diverse needs of clinics and hospitals.



LGT-3500LP

Upright Design



LGT-3500GP

Desktop Design



DUKON ELECTROTHERAPY CONT..

Based on neurophysiology,magnetobiology, biophysics, and clinical brain disease treatment, a special treatment device outputs specific regular transcranial magnetic stimulation. The physiological effects of transcranial magnetic stimulation on brain tissue are mainly manifested in the following aspects:

- Increase the elasticity of cerebral blood vessels.
- Improve blood circulation in the lesion area. mprove the metabolic environment of brain tissue.
- Guide the patient's brain magnetic function towards normalization and orderliness



MAIN TECHNICAL PARAMETERS			
rTMS Intensity	Magnetic Field Frequency	Erebellar fastigial nucleus stimulation Maximum Output Current	Limb electrical stimulation
III LEVEL 3-18mT	50Hz±1Hz	< 20mA	60mA±30%



PC version for double use



Touch screen for single use/dual channels



Touch screen for single use (Optional trolley)



Home version

DUKON ELECTROTHERAPY

REPETITIVE TRANSCRANIAL MAGNETIC STIMULATOR

The repetitive transcranial magnetic stimulator uses the principle of electromagnetic induction to discharge rapidly from the energy storage capacitor to the stimulation coil. The magnetic field generated by the coil current induces the biological current in the biological tissue of the stimulation site through the skull or other tissues, which can stimulate the brain or peripheral nerves harmlessly and painlessly, causing the excited activity of nerve cells. It can be used not only for the diagnosis of nerve function, but also for the intervention and treatment of nerve function regulation. It is a non-invasive and accurate brain science and neuroregulation technology.




8-Shaped coil
(Standard)

0-Shaped coil
(Standard)

Therapy cap
(Optional)

MAIN TECHNICAL PARAMETERS

Magnetic Induction Intensity	Output Frequency	Pulse Width	Output Time	Low Magnetic Intensity (Optional)
0-2T($\pm 20\%$) divided into 100 levels	0-30Hz($\pm 10\%$) in steps of 1Hz	350 μ s ($\pm 10\%$)	0-120mins	3-18mT ($\pm 10\%$)



★ PRODUCT ADVANTAGES

- Integrated computer touch screen with an intelligent and simple human-machine interface for convenient operation
- Consistent stimulation output(intensity, amplitude, etc.) without decay during continuous stimulation to ensure stable and reliable stimulation quality
- Multiple temperature monitoring,unique cold processing system, intelligent liquid cooling, and dual air cooling functions for safe and long-duration continuous operation. Ensures patient safety
- Pure waveform, free from harmonics and interference, resulting in higher treatment precision
- Fuselage design emergency stop button, users can stop treatment at any time during use

DUCON

SWALLOWING NEUROMUSCULAR ELECTRICAL STIMULATOR

Swallowing neuromuscular electrical stimulation devices are used for the treatment and training of swallowing dysfunction. They output specific low-frequency pulse currents to electrically stimulate the nerves and muscles in the laryngeal and neck area, exciting the nerves and swallowing muscle groups, promoting the recovery of damaged nerves, strengthening the movement of the swallowing muscles, alleviating neuronal paralysis, and promoting the reconstruction and restoration of the swallowing reflex arc function, thereby improving swallowing and language abilities.



★ MORE MODELS



**Dual Channel
Computer Version**
(Includes hand controller)



**Dual Channel
Touch Screen Version**
(Includes hand controller)



**Single Channel
Touch Screen Version**
(Includes hand controller)



**Single Channel
Desktop Version**
(Includes hand controller)



Home Portable
(optional hand controller)

LOWER LIMB FUNCTIONAL ELECTRICAL STIMULATOR

Lower limb functional electrical stimulator use electrical stimulation to mimic the central nervous system, controlling targeted nerves and muscles. These devices are worn by patients, and their built-in motion sensors detect the patient's gait, assisting them in completing walking actions such as knee flexion, foot lifting, stepping, and landing. The devices can increase joint mobility, alleviate lower limb spasticity, improve gait, and reduce the risk of falls.



★ PRODUCT ADVANTAGES

- It integrates passive electrical stimulation, electromyographic feedback stimulation, and multimedia biofeedback training
- Three-dimensional gait analysis training provides real-time, intuitive feedback of rehabilitation indicators to patients, which is highly effective in improving gait



Main
unit



Foot sole
transmitter



Remote
control


ELECTRICAL DIAPHRAGM PACEMAKER DUKON

LOW-FREQUENCY EXTERNAL DIAPHRAGM PACEMAKER

By using Skin electrodes to apply low-frequency pulsed electrical stimulation to the phrenic nerve, the diaphragm contracts and relaxes rhythmically, increasing diaphragmatic excursion, thereby increasing ventilation, promoting the expulsion of CO₂ from lungs, and gradually restoring the patient's diaphragmatic function. It's an important equipment for critical care rehabilitation and cardiopulmonary.



Electrode lead wires



Electrode pads

MAIN TECHNICAL PARAMETERS			
Pulse Repetition Frequency	Pacing Times	Treatment Time	Output Amplitude
30Hz, 35Hz, 40Hz, 45Hz, 50Hz	Adjustable from 5 to 15	Adjustable from 5 to 90 minutes	6mA (adjustable in 60 levels)

ELECTRICAL STIMULATOR DUKON

MEDICATED PERMEATION THERAPY DEVICE

Medicated permeation therapy device is a medical equipment that integrates the wisdom of traditional Chinese medicine and the essence of modern science and technology. Combining medium frequency electrotherapy and drug penetration technology, it realizes the synergistic effect of medicine and physical therapy. The modulated wave is sent to the human body with the strong penetrating power of the carrier if current to achieve the effect of deep treatment, and the drug can penetrate into the lesion site more quickly by using the pushing or heating effect of the IF current to achieve accurate and efficient treatment. The equipment is easy to operate, safe and non-invasive.




Single channel



Dual channel

★ PRODUCT ADVANTAGES

The product is composed of three parts: medium frequency electrotherapy, warm electrode and three-dimensional dynamic interferometer, which can realize multiple functions such as medium frequency physiotherapy, medium frequency drug introduction, interferometer drug introduction and interferometer electric physiotherapy. It is effective in relieving inflammation and analgesia, improving blood circulation, promoting the regression of inflammation and edema, accelerating fracture healing, preventing and delaying muscle atrophy.



ELECTRICAL STIMULATOR DUKON

INTERFERENTIAL CURRENT THERAPY DEVICE



The device is three different frequency equal amplitude alternating current cross input human body at the same time, produce a constantly changing interference current; Can stimulate the body from multiple angles and directions;

Thus, the dynamic stimulation effect of multiple parts, different directions, angles and shapes can be obtained. The three dimensional interference electric field was generated in the body during treatment, and the analgesic and analgesic effects were obvious.

★ PRODUCT ADVANTAGES

- There are a variety of different frequencies and waveforms of current, can be flexibly selected according to the patient's condition, the patient is not easy to produce tolerance
- With fault alarm and other functions, the operation is simple and safe



Single screen Single channel



Two-screen Two-channel



Portable Single channel

FUMIGATION THERAPY DUKON



Fumigation therapy is one of the commonly used external treatment methods in traditional Chinese medicine and is an important part of Chinese medical science. The fumigation therapy device mainly achieves therapeutic purposes through the dual effects of physical warmth and absorption of traditional Chinese medicinal herbs. By using heat to apply the medicine to the steamed area, it can expand local blood vessels, promote blood circulation, warm and unblock the meridians, detoxify and kill bacteria, relieve itching, reduce swelling and pain, and ultimately achieve the goals of treating diseases, preventing diseases, health care, and beautification.

★ PRODUCT ADVANTAGES

- Real-time temperature control design, effectively prevent skin burns
- Universal arm design, fumigation Angle 360° adjustable, convenient and flexible
- Anti-condensation water dripping design to avoid skin discomfort



Available in single channel and dual channel specifications

CHEST VIBRATOR DUKON

CONSTANT TEMPERATURE WAX THERAPY DEVICE



- Automatic control, melting wax, washing wax, filtration, automatic wax and wax cake production, no need to manually connect wax or scoop wax.
- The combination of ultraviolet light and high temperature disinfection can effectively kill the bacteria on the surface and inside of the paraffin to prevent cross infection and ensure safer use.
- Paraffin cleaning: can wash the reused paraffin for separation, automatic disinfection, cleaning, precipitation, multistage filtration, to ensure that the paraffin is clean and health can be reused for a longer time, and completely solve the problem of paraffin cleaning in the traditional wax therapy machine.
- Special air duct circulation system design of wax cake box ensures that the temperature difference at each point of wax cake box is not more than $\pm 0.5^{\circ}\text{C}$ when the wax cake box is fully loaded, to ensure the coexistence of no sandwich and no wax liquid inside wax cake, to ensure the softness, adhesion and plasticity of wax cake, and to ensure safe and effective use.
- 10" super large screen display interface, independent touch algorithm design, clearer display, more smooth operation.
- The automatic model has 12/18 discs to choose from, with a single disc waxing system, and the wax cake thickness can be adjusted at three levels.

WET HOT COMPRESS THERAPY DEVICE



Wet hot compress therapy device is the first domestic use of glass fiber and pure materials, through high-quality cotton and linen canvas packaging, heating in the water tank and absorb a lot of heat and water, the treatment of slow heat dissipation and emit water vapor, the treatment site play a role in wet hot compress. The hot bag can be used repeatedly until the material fails. The device is simple to operate, and the treatment site without "pollution" and unclear feeling, with pure natural non-toxic side effects, ultra-long heat preservation, repeated durability and other characteristics, is an economical and practical new thermal therapy product.

MULTI-FREQUENCY VIBRATION SPUTUM DISPENSER



According to the principle of clinical chest physical therapy, the multi-frequency vibration expectorator produces a therapeutic force with a specific direction of periodic changes on the surface of the patient's body, specifically for patients with lung ventilation disorders. Simulating the artificial back tapping, there are percussion, tremor and pushing mode of work, fully realize the purpose of drainage for patients and postures.

LASER THERAPY HAOBRO

★ PRODUCT PRINCIPLE

Magnetic therapy: Strong magnetic fields can adjust the biological magnetic field of the body, generate induced microcurrents, change the permeability of cell membranes, change the activity of certain enzymes, expand blood vessels, and accelerate blood flow, thereby achieving auxiliary therapeutic effects such as pain relief and swelling.

Laser: It is mainly used to produce photothermal, photochemical and biological stimulation effects with human tissues to achieve the purpose of treatment or auxiliary treatment. After biological tissue absorbs the energy of low energy laser, its energy state changes, accompanied by few thermal effect and photochemical effect to stimulate the tissue. By improving blood circulation and promoting the dissipation and absorption of inflammation, it is beneficial to tissue regeneration and accelerates wound healing.



★ FEATURES

- Extra large touch screen**
One-button shuttle design makes operation easy and convenient
- 3 working modes**
Single trigger, continuous mode, automatic mode
- 6 Tesla high intensity magnetic field**
Strong penetration and deep treatment depth
- Non-invasive treatment**
Painless, non-invasive, good experience
- Using 3 laser sources**
655nm wavelength low energy laser
- Liquid cooling system**
Long working time and uninterrupted treatment
- Built-in commonly used clinical treatment prescription**
Also supports customized treatment prescription
- Multiple protection**

TECHNICAL SPECIFICATIONS

BC-MT01	BC-MT02	BC-MT03	BC-MT04	BC-MT05
Cart-based type	Cart-based type	Cart-based type	Cart-based type	Cart-based type
Single channel	Dual channel	Single channel	Dual channel	Dual channel
Treatment head 1	Treatment head 1 +treatment head 2/3	Treatment head 1	Treatment head 1 +treatment head 2/3	Treatment head 3
/	/	1 workstation	1 workstation	1 workstation



OSTEOPOROSIS – HAOBRO

★ PRODUCT PRINCIPLE

When a specific low-frequency pulsed magnetic field acts on the human body, it will produce a magnetoelectric effect. After acting on the bone cells, it improves the metabolism and reconstruction of the bone, promotes the activity of the osteoblast, and delays the loss of bone mass. At the same time, magnetic field analgesia can relieve pain caused by osteoporosis.



★ FEATURES

1. Double ring magnetic therapy device
2. Planar magnetic therapy device (7-section crawler type)
3. Carbon fiber infrared heating device
4. Narrow gap design, beautiful and low failure rate
5. The treatment ring is driven and positioned by a precision screw
6. Six treatment modes

TECHNICAL SPECIFICATIONS

HB320 Double ring stereo type

1 double ring stereo treatment bed

Bed surface heating device

HB330 Double ring stereo type

2 double ring stereo treatment beds

Bed surface heating device



BONE THERAPY – HAOBRO

★ PRODUCT PRINCIPLE

The electric current is introduced into the human body to stimulate the growth of bone cells; the pulsed magnetic field acts on bones, muscles, soft tissues, etc. to produce a series of magnetic biological effects; and through the electrostatic field, the piezoelectric effect of the bone is formed, which promotes the deposition of calcium salts and accelerates calcification, so as to treat bone tissue damage.

★ FEATURES

1. Comprehensive treatment of various physical factors, magnetic therapy, low-frequency electrotherapy, electrostatic therapy
2. Magnetic field strength up to 50mT
3. Paramagnetic, focusing and other magnetic therapy modes
4. Multiple channels can be used at the same time
5. The magnetic coupling disc can be adjusted arbitrarily, which can meet the treatment needs of various parts



TECHNICAL SPECIFICATIONS

DHG-I	DHG-III
Cart-based	Cart-based
Touch screen	Touch screen
One-key shuttle	One-key shuttle
2 channels (magnetic therapy)	1 channel (magnetic therapy)
4 channels (low frequency electrotherapy)	2 channels (low frequency electrotherapy)
2 channels (electrostatic therapy)	/

HAOBRO ULTRASONIC DEVICE

★ FEATURES

1. 8-inch LCD screen design, real-time display of various treatment data
2. One-key shuttle button for easy operation
3. One-key switch between pulse and continuous mode, and the treatment intensity can be adjusted in 10 levels
4. Four channels can be used at the same time
5. Combination of fixed treatment head and mobile treatment head



TECHNICAL SPECIFICATIONS			
HB810A	HB810B	HB810D	HB820D
Single channel	Dual channel	Dual channel	4-channel
Portable	Portable	Cart-based	Cart-based
/	One 1MHz fixed treatment head	One 1MHz fixed treatment head	Two 1MHz fixed treatment heads
One 1MHz mobile treatment head	One 3MHz mobile treatment head	One 3MHz mobile treatment head	A 1MHz and a 3MHz mobile treatment head respectively



ULTRASONIC PHYSIO THERAPY DEVICE

★ FEATURES

1. Color touch screen+knob operation
2. Multiple ultrasound treatment heads available for selection
3. Ultrasonic treatment head with over-temperature protection function
4. Built in multiple clinical treatment plans
5. Meeting different clinical needs
6. Two ultrasonic operating frequencies 1MHz, 3MHz dual frequency output, single treatment head can realize dual frequency free switching
7. Customized parameter collection function
It can collect 1000 custom solutions
8. Optional rechargeable battery
It can be used normally without AC power supply
9. With an overview interface to view the operating status of each channel in real time
10. Quick guidance plan
Generate treatment parameters based on treatment depth, tissue status, and treatment area



TECHNICAL SPECIFICATIONS

HB-UL1	HB-UL2	HB-UL3	HB-UL4
1 pluggable ultrasonic channel	2 pluggable ultrasonic channels	1 adsorption ultrasonic channel	2 adsorption ultrasonic channels
Standard with A28 type	Standard with A28 type, A13 type	Standard with A28V type	Standard with A28V type
Optional A13, A25,	Optional A25, A28V,	/	/



RADIAL SHOCKER THERAPY HAOBRO

★ PRODUCT PRINCIPLE

The compressed air generated by the program control of the pulse acoustic wave into a highly accurate ballistic shock wave, the use of equipment unique frequency, through the positioning and movement of the Shockwave Transmitter can produce cavitation effect, stress effect, analgesic effect, metabolic activation effect on the human body.

★ FEATURES

1. 10.1-inch color touch screen
2. Built-in and customized protocols
3. 4 pain assessment function modes with evaluation results
Dynamic VAS, static VAS, sleep VAS, facial expression measurement chart
4. Patient data management
Can input patient information, view, edit, delete, search, add, and load
5. Protective function
Handpiece detachment alarm
Insufficient air pressure alarm
6. Step intensity and step frequency can be set
7. The service life of the transmitters and bullet is up to 2 million shots
8. Two types of transmitters

TECHNICAL SPECIFICATIONS

HB100M	HB101M	HB102M
Single channel	Single channel	Dual channel
Portable	Cart-Based	Cart-Based
Shockwave Handpiece x1 +Vibration Handpiece x1		Shockwave Handpiece x2 +Vibration Handpiece x1



ELECTRO MAGNETIC HAOBRO

Using computer information simulation technology to synthesize simulated biological current and biological magnetic field, through the electrode pads attached to the mastoid parts on both sides and special magnetic therapy caps, it can non-invasively penetrate the skull barrier from outside the skull to stimulate the parietal nucleus area and improve the elasticity of cerebrovascular, increase local blood flow in the brain, promote brain microcirculation, accelerate the metabolism of brain cells, and accelerate the recovery of brain functions.



★ FEATURES

1. Four powers in one
2. Dual screen display, equipped with a shuttle button, easy to operate
3. Built-in 11 kinds of treatment prescriptions
4. Four unique treatment frequencies
23.81Hz, 15.87Hz, 15.87Hz, 11.90Hz
5. Four treatment modes

TECHNICAL SPECIFICATIONS

HB510B	HB520B	HB510D	HB520D
Single channel	Dual channel	Single channel	Dual channel
Portable	Portable	Cart-based	Cart-based
1-channel main electrode treatment	2-channel main electrode treatment	1-channel main electrode treatment	2-channel main electrode treatment
2-channel auxiliary electrode treatment	4-channel auxiliary electrode treatment	2-channel auxiliary electrode treatment	4-channel auxiliary electrode treatment
1-channel rTMS treatment	2-channel rTMS treatment	1-channel rTMS treatment	2-channel rTMS treatment



ELECTROTHERAPY NERVE STIMULATION HAOBRO

Nerve and muscle electrical stimulators output specific low-frequency pulse currents, which can directly excite nerves and muscle tissues, cause rhythmic muscle contractions, and improve muscle strength.

★ FEATURES

1. Color touch screen & one-key shuttle design, easy to operate
2. Electroacupuncture treatment
3. Negative pressure suction therapy
4. Electrode therapy



TECHNICAL SPECIFICATIONS			
HB-SJ1	HB-SJ3	HB-SJ5	HB-SJ6
Single channel dual output	Dual channels and four outputs	Three channels and six outputs	Eight channels and eight outputs
Portable	Portable	Cart-based	Cart-based
1-channel electrode + 1-channel electric needle output	2-channel electrodes + 2-channel electric needle outputs	2-channel electrodes + 3-channel electric needles + 1-channel negative-pressure output	8-channel electrodes



SPASMODIC LOW STIMULATOR HAOBRO

Two sets of specific low-frequency pulse currents are alternately output – one set is used to stimulate the spastic muscles, and the other set of currents stimulate the antagonist muscles. Through the mutual inhibition of the two sets of currents, the spastic muscles are relaxed, the muscle tension is reduced, and the muscle strength of the antagonist muscles are increased, thereby the functions of the limbs is improved.

★ FEATURES

1. Color touch screen & one-key shuttle design
Real-time display, easy to operate
2. Two treatment electrodes
Repetitive treatment electrodes; disposable treatment electrodes
3. Negative pressure suction therapy
4. Built-in 10 kinds of treatment prescriptions



TECHNICAL SPECIFICATIONS

HB-JL1	HB-JL3	HB-JL5
Single channel and dual outputs	Dual channels and four outputs	Three channels and six outputs
Portable	Portable	Cart-based
2-channel electrode output	4-channel electrode output	4-channel electrode + 2-channel negative-voltage outputs



MYOELECTRICAL STIMULATOR HAOBRO

By collection of the patient's surface myoelectricity signal, use of electronic information technology to amplify it, and conversion of it into visual and auditory feedback to the patient and medical staff, so that the patient can understand their own physiological changes, and through conscious changes, the patient's physiological changes are directed towards a certain target direction for improvement to carry out physical and mental illness and physical rehabilitation.

TECHNICAL SPECIFICATIONS

HB120B	HB120D
Dual channel	Dual channel
Portable	Cart-based



★ FEATURES

1. 8-inch LCD touch screen display
2. Electrode therapy
3. Built-in multiple treatment modes
4. Built-in 61 treatment prescriptions

★ CLINICAL ADVANTAGE

- It collects and analyzes biological myoelectricity through software without subjective judgment of medical staff;
- After obtaining the results through software analysis, it uses low-frequency electrical stimulation of nerves and muscles to restore the function of the stimulated nerves and muscles;
- It is of high efficiency and can reduce the work intensity of medical staff;
- It allow patients to clearly understand the status and changes of their physiological activities and increase their confidence.

SHORT WAVE THERAPY HAOBRO

The method of treating diseases using high-frequency electromagnetic fields generated by shortwave current is called shortwave therapy. The short wave therapy device can generate a clinical electromagnetic field with a frequency of 27.12MHz, which acts on the human body through the application device, generating different levels of heat at the treatment site to treat related diseases.

TECHNICAL SPECIFICATIONS

HB-SW200C	HB-SW200L	HB-SW400
Vertical	Vertical	Vertical
Capacitive applicator	Inductive applicator	Capacitive applicator + Inductive applicator



★ FEATURES

1. **Portable operation**
LCD touch screen + one key shuttle
2. **Two kinds of output**
Continuous output, pulse output
3. **Automatic tuning**
Bring the treatment to an appropriate output state
4. **Diversification**
Applicator selection
5. **High stability**
High energy conversion rate and low radiation interference

ELECTRIC DEVICES HAOBRO

Apply a certain frequency of current to the body, to cause comfortable muscle contractions, to train muscles, to inhibit pain fiber conduction, to release painkillers, to improve blood circulation, to soften the scar, to release the conglutination.

★ FEATURES

- 1. 7 / 10.1" color touch screen**
With aesthetic Shape, easy to operate
- 2. Multiple interference modes**
With normal, dynamic and other modes, apply for different conditions, meet different clinical needs
- 3. Built-in treatment prescription**
Provide professional treatment pro-grams,

Provide suggestions and custom-ized prescription functions, convenient for clinical operation

- 4. Multiple security protection**
- 5. Electrode sheet heating function**
Electrode sheet heating function, human-ized design, improve human comfort



TECHNICAL SPECIFICATIONS

HB-ZP10	HB-ZP30	HB-ZP40	HB-ZP50
Portable	Portable	Cabinet Type	Cabinet Type
2 channels	4 channels	4channels	8channels
IF+ interference	IF+ interference	IF+ interference +iontophoresis	IF+ interference +iontophoresis



DRUG INTRODUCTION THERAPY HAOBRO

★ FEATURES

- 1. 5-inch color LCD touch screen**
- 2. In hand-held portable design, it can be freely combined into a desktop model**
- 3. With counting function, it can automatic identify the counting card and read the number of uses**
- 4. Hyperthermia function: Temperature adjustment in the adjustment range of 38°C ~ 44°C**
- 5. Three treatment modes, multiple prescription options, and five output waveforms**



TECHNICAL SPECIFICATIONS

HB-LZ21SW	HB-LZ22SW	HB-LZ21S
2 (1 pulse + 1 ultrasound)	2 (2 pulses)	2 (1 pulse + 1 ultrasound)

INTERFERENCIAL ELECTRICAL DEVICER HAOBRO



Cross input of two groups or more of intermediate frequency currents of different frequencies into the human body, interfering with each other in the body to generate difference frequency currents, which can relieve pain, excite neuromuscular tissue, promote local blood and lymphatic circulation, loosen adhesion, and improve smooth muscle strength, etc. effect.

★ FEATURES

1. 7-inch true color touch screen display
Good-looking appearance and convenient operation
2. Dual-channel output mode
Multiplex negative pressure adsorption therapy, Multiplex electrode therapy
3. Five kinds of carrier frequency
4. Ten interference modes, suitable for different diseases and meeting different clinical needs
5. Six kinds of dynamic rhythms are adjustable: 0~15s, tolerance $\pm 3s$
6. Top plate heating function
Humanized design increases comfort

TECHNICAL SPECIFICATIONS

HB-ZP2	HB-ZP4	HB-ZP6
Single-channel countertop	Single-channel cabinet type	Dual-channel cabinet type
3-channel electrodes + 3-channel negative pressure outputs	3-channel electrodes + 3-channel negative pressure outputs	6-channel electrode + 6-channel negative pressure outputs
Stereo dynamic interference	Stereo dynamic interference	Stereo dynamic interference
Top plate heating	Top plate heating	Top plate heating

MICROWAVE THERAPY HAOBRO



The 2450MHz high-frequency electromagnetic field acts on the affected area to cause high-frequency oscillations of ions, water molecules, and dipoles in tissue cells to generate heat, which can promote blood circulation, improve metabolism, accelerate tissue repair, reduce inflammation, reduce swelling, relieve pain, and relieve muscle spasm, etc.

TECHNICAL SPECIFICATIONS

AMT-BI	AMT-C
Cart-based	Cart-based
Microcomputer control, single microwave source	Microcomputer control, double microwave source



★ FEATURES

1. The operation is simple and can be mastered by general medical staff
2. The operation time is short, which is usually only a few seconds
3. The patient has no pain or bleeding during the operation
4. There is no sputtering to medical staff and no smell and fog of tissues burnt
5. The treatment effect is good, the time required for healing is short, and there is no scar

THERMAL TENS THERAPY HAOBRO

Thermotherapy is the use of water to transfer heat energy to the body, producing biological effects in the body to achieve the goal of preventing and treating diseases. Low frequency pulse current stimulates nerve or muscle contraction to achieve therapeutic effects, causing muscle recruitment activity, enhancing muscle strength, improving local blood circulation, promoting the recovery of central motor control function, and reconstructing normal motor modes.



★ FEATURES

1. 7-inch color touch screen with rotary button design for efficient operation
2. Output electrode with one positive and two negative, with left and right balance adjusted as needed
3. Warm function with six options
4. Multiple treatment modes available for clinical selection
5. Built in negative pressure suction

TECHNICAL SPECIFICATIONS

SL-RZB-A	SL-RZB-C
Vertical type	Portable type
2 channels	1 channel
Warm electric stimulation+negative pressure suction	Warm electric stimulation+negative pressure suction



STEAM THERAPY HAOBRO

★ FEATURES

1. Five utility model patents
2. Smart countdown function high stability and safety
3. Water and electricity isolation design



TECHNICAL SPECIFICATIONS

HB1000	HB4000
Single nozzle	Double nozzle

DEEP MUSCLE STIMULATOR HAOBRO

✦ PRODUCT PRINCIPLE

The continuous vibration of the instrument acts on deep muscles and soft tissues, effectively alleviating tissue adhesion and excessive excitement, quickly decomposing inflammatory products, promoting metabolism, and improving the tissue environment; at the same time, the receptors are constantly affected by vibration stimulation, thereby pain transmission is inhibited.

✦ FEATURES

1. All-metal treatment device of strong power, imported from Europe
2. In portable handbag, the device is small and light
3. The treatment probe is made of titanium alloy and has good mechanical properties
4. Innovative design, with three specifications of treatment probes, meeting various clinical treatment needs of human myofascial chain
5. Multi-parameter adjustable (HB-DMS200): vibration frequency 15-60Hz, treatment time 0-30min
6. The amplitude of the treatment head is 6mm and the highest can hit 3,600 times per minute



TECHNICAL SPECIFICATIONS

HB-DMS100 (fixed frequency)

Frequency: 60Hz Amplitude: 6mm

HB-DMS200 (variable frequency)

Time: 0-30min Frequency: 15-60Hz Amplitude: 6mm

Titanium alloy probe: 20mm, 25mm, 35mm

Vibration meter size: 270mm x 150mm x 54mm

CHEST WALL EXPECTORATOR HAOBRO



✦ PRODUCT PRINCIPLE

The host generates high-frequency air pulses through rapid inflation and deflation. The air pulses inflate the airbag vest through the catheter to generate vibrations, thereby compressing the chest cavity at a certain frequency, causing the airflow to quickly pass through the lung lobes and respiratory tracts, and promoting mucus and metabolites on the mucosal surface to get loose, be liquefied, fall off, and move to the mouth. After reaching the mouth, phlegm is discharged through normal cough or negative pressure suction.

✦ FEATURES

1. High-definition touch screen, easy to operate
2. Dedicated modes for adults and children respectively
3. One-button emergency stop, stable and reliable
4. Two functions in one: sputum expectoration and atomization
5. Especially made vests of multiple models

EXPECTOLATOR VIBRATION CHEST HAOBRO

★ PRODUCT PRINCIPLE

The percussion head produces vertical and horizontal therapeutic forces on the surface of the human body. The percussion and tremor produced by the vertical therapeutic force can promote the relaxation and liquefaction of mucus and metabolites on the mucosal surface of the respiratory tract, and the directional push generated by the horizontal force helps the liquefied mucus to be excreted from the body according to the selected direction.

★ FEATURES

1. 11 different specifications of treatment heads
2. 10 treatment modes (for adults and children)
3. Two-in-one: sputum expectoration and atomization
4. Super large touch screen, simple and easy to operate interface display
5. With hand-held body position drainage operation, easier to use
6. Satisfying the need for expectoration in multiple directions vertically and horizontally



TECHNICAL SPECIFICATIONS

HB-KPT4	HB-KPT6
Dual channel	Dual channel
Cart-based	Cart-based
Sputum drainage	Expectoration + atomization
Adults and children	Adults and children
220VA	3000VA

LOWER LIMBS VIBRATION TRAINING HAOBRO

★ FEATURES



1. Console operation

Connect the power cord and turn on the power switch. The control console is used to control the treatment parameters and manually adjust the treatment bed placement angle. The treatment parameters are intensity, frequency, preset parameter mode and treatment time. Under the control of the console, the vibration machine can achieve various different frequencies of regular vibration to meet the training needs of different parts of the human body.

2. Forearm pallet height/rotation operation

The height and angle of the forearm bracket can be adjusted according to the actual needs of the user to meet the use of different age levels of users.

3. Emergency stop switch

In the training process, if you need the product emergency stop, you can directly press the emergency stop button, as shown in the figure. When the emergency stop button is pressed, the red light is on, indicating that the product is stopped. Press the emergency stop button again and the product will resume working.

★ PRODUCT PRINCIPLE

The lower limb vibration training system consists of a base, a machine base, a treatment bed, a vibrating machine, and a control console.

INFRARED THERAPY HAOBRO

PRODUCT PRINCIPLE

Infrared therapy

When infrared rays act on the body, the skin and subcutaneous tissues convert the absorbed infrared energy into heat, which can cause vasodilation, improve local blood circulation, promote metabolism, reduce muscle tension, and reduce the excitability of sensory nerves, thereby achieving therapeutic effects such as analgesia, detumescence, and relieving muscle spasms

Polarized light therapy

Polarized light therapy, based on traditional meridian theory, directly irradiates acupoints with small light points; Guided by modern theory, direct irradiation is applied to local lesions and ganglia, nerve trunks, nerve plexus and other parts, and the combination of the two theories can achieve the therapeutic effect of light acupuncture and light moxibustion.



FEATURES

1. True color touch screen display
2. Combining points and surfaces for a wider range of treatment options
3. Three channel independent control
4. Multi directional and angle adjustable bracket

TECHNICAL SPECIFICATIONS

SL-PZG-A	SL-PZG-B	SL-PZG-C
LCD display	LCD display	LCD display
1 divergent infrared radiation head +2 clustered polarized light radiation head	2 Clustered Polarized Light radiation head	1 divergent infrared radiation head +1 cluster polarized light radiation head



INFRARED THERMAL THERAPY HAOBRO



★ FEATURES

1. Large true-color touch screen

The operation menu is simple and clear

2. Dual channel independent control

It can treat two patients at the same time

3. Powerful power and deep treatment

High power light source: optical density > 235mW/cm² at 20cm

4. Four safety protection measures

Power failure protection, anti-tipping protection, overheat alarm, safety distance setting

5. Electric lifting system

The treatment angle can be adjusted and easily positioned according to the patient's body position

6. Patient information management and memory function

With patient information management, treatment plan and memory function

★ PRODUCT PRINCIPLE

The infrared thermal radiation therapy device mainly outputs infrared therapeutic irradiance for the treatment of pain and inflammation, which can improve blood circulation, promote tissue repair and regeneration, eliminate swelling, and accelerate wound healing.



HB-HW750	HB-HW750A	HB-HW750B
Single light source	Single light source	Dual light source
Manual lifting	Electric lifting	Electric lifting

SMART MASSAGE THERAPY HAOBRO



★ FEATURES

1. A new method of pain management.

Muscle massage, chiropractic, deep muscle stimulation and other functions are available. The combined therapy has good curative effect and high efficiency

2. 7-inch color touch screen

7-inch color touch screen design, easy to operate

3. The vibration speed of the massager can be adjusted in 9 gears, and the hitting force of the chiropractic gun

can be adjusted in 10 gears

Meet the needs of various clinical treatment intensity output

4. A variety of massage head and chiropractic head optional

Increase departmental revenue and meet different clinical needs

5. Multiple treatment modalities

Manual, automatic, cervical spine, thoracic spine, lumbar spine, limb joints and other modes can be selected, humanized design, easy to clinical

★ PRODUCT PRINCIPLE

AM series Smart Massage Chiropractic device uses horizontal and vertical rhythmic vibration and adjustable pulsed high-energy impact force. The two therapies are combined to work on the deep muscle and skeletal system, thereby improving the level of microcirculation metabolism and bone structure. It has a good auxiliary effect on muscle spasm, soft tissue pain and spinal correction respectively.

RHYTHM THERAPY DEVICE- FIZYOMED



★ FIELDS OF APPLICATION


- Tendon damages
- Muscle spasms
- Spinal and Cervical disc hernia and pains
- Muscle tearing
- Cerebral palsied
- Cellulite
- Capillary vessel system diseases (varicosity)
- Sclerosis
- Weight Loss treatment
- Tendon injuries
- Osteoporosis pains
- Lymph and edema treatments

Fizyoterra 1000 rhythm therapy device applies mechanical vibrations identical to natural physiological vibrations in human body by means of the particular therapy headpiece and it stimulates the nervous system and muscles rhythmically. As a result of these vibrations applied between 8 -12 Hz. relaxing effect is experienced by the patients. Consequently, metabolism is returned to its normal state in a short time and recovery and renovation process in the tissues restarted at cellular level.

TECHNICAL SPECIFICATIONS:

Güç Girişi / Power Supply	220 - 240V 50/60 Hz
Ölçüleri / Dimensions	200x250x210 mm
Ağırlık / Weight	5 kg
Tedavi Zamanı / Therapy Duration	1 - 99 dk / min
Hız / Speed (Max)	%10 - %80
Sınıf Beyanı / Class Declaration	class I
Sertifikalar / Certificates	EN ISO 13485 : 2016
	CE (Declaration Of Conformity)

INFRARED UNIT-FIZYOMED



★ FEATURES

- 2 treatment lamps
- Front end of the unit is protective.
- Adjustable height
- Single use or dual use of the treatment lamp is available
- Treatment time adjustable
- Easy mobility
- Heat resistance electrical cabling
- Anti rust electro static paint on all metal / steel parts

TECHNICAL SPECIFICATIONS

Boy / Legth	130-200 cm arası ayarlı
Eni (Tekerlekli Ayak Kısmı Dahil) / With Of The Leg With Castors	46 cm
Ampul Şiddeti / Bulb Intensity	2 x 250 Watt
Tedavi Zamanı / Therapy Duration	0-60 dk/min
Yaklaşık Ağırlığı / Approx. Weight	8 kg
Sınıf Beyanı / Class Declaration	class I
Sertifikalar / Certificates	EN ISO 13485:2016
CE (Declaration Of Conformity)	min. 80
TPE Tescil no: 2008 02295	125/250
Standart Aksesuarları / Standard Accesories	2 Adet 250 Watt Infraruj Lambası

INVACMED-MED INMED



★ FEATURES

- 2 channels.
- Connections for pneumatic – electric cables with vacuum electrodes.
- Continuous or pulsed vacuum mode.
- Fast, practical fixing of electrodes in hard-to-reach places of the body.
- Perfect electrode adherence to the body.
- Cooperation with electrotherapy units.
- Possibility to applicate the currents through vacuum electrodes (with integrated metal electrodes).
- Connections for conventional electrodes for electrotherapy.
- Electronic control and protection against too high vacuum force.
- Indication of water tank filling.
- Power saving operation.
- Display located at an angle making it easier to operate the device

TECHNICAL DATA:

Pulse frequency (pulsed mode) [pulses/min.]:	15-90
Vacuum [mbar]:	60-400
Number of electrodes:	4
Power supply [V/Hz]:	230/50-60
Current consumption [A]:	0,1
Dimensions [cm]:	33,5 x 13,5x 32,5
Weight [kg]:	3,4
Dehydrator reservoir capacity [ml]:	min. 80
Fuse [mAT/V]:	125/250
Optical alarm level "H ₂ O"	no more than 50 ml



MAGNETRONIC MF-2,MF-12-MED INMED



★ FEATURES

- 2 magnetotherapy channels.
- Possibility of performing two treatments at the same time.
- Colour touch screen display (4.3").
- User-friendly touch screen and button controls.
- Ready to use preset treatment programs for popular treatments.
- User-defined programs can be stored by the therapist.
- Individual adjustment of treatment settings.
- Fan control feature to minimize noise and power consumption.
- Treatment duration and number of performed treatments counter.
- Can be used as portable device.
- Magnetic field shapes: sine, rectangle and triangle, each of them unipolar or bipolar,
- MX1 - consecutive shape change at constant frequency, or
- MX2 - consecutive shape change at variable frequency

TECHNICAL DATA:

Magnetic field frequency range [Hz]:	1-150
Magnetic field intensity [mT]:	0-8 (depends on applicator used)
Pulse / break duration [s]:	0.5-8
Work mode:	continuous or interrupted
Treatment timer [min]:	1-99
Power supply:	~230V 10%, 50Hz, 70VA
Electrical protection class:	I, type B
Dimensions (L x W x H) [cm]:	33,5 x 27 x 12,5
Weight (of control unit) [kg]:	2,5



★ FEATURES

- 2 independent treatments. The device features an innovative control system – each applicator has an independent treatment timer and can be independently initiated. As a result, Magnetronic MF-12 can perform 2 independent magnetotherapy treatments in different time, with the same treatment settings.
- Connections for two applicators

TECHNICAL DATA:

Magnetic field frequency range [Hz]:	1-100
Magnetic field intensity [mT]:	0-20
Pulse/break duration [s]:	0,5-8
Power supply [V/Hz/W]:	230/50/200
Weight of control unit [kg]:	4,4
Dimensions of control unit (L x W x H) [cm]:	33,5 x 27 x 12,5
Shapes of magnetic field modulation:	sine, rectangle, triangle – unipolar and bipolar

MAGNETRONIC MF24-MED





★ FEATURES

- 5 simultaneous treatments
- The device features an innovative control system – each applicator has an independent treatment timer. As a result, Magnetronic MF-24 can perform 4 simultaneous independently initiated magnetotherapy treatments (2 treatments per channel) and 1 laser therapy treatment. In total Magnetronic MF-24 features 5 treatment timers.
- 2 independent magnetotherapy channels with separate adjustments.
- Each channel can connect 2 applicators.
- Independent laser therapy channel for treatments with a laser probe.

TECHNICAL DATA:

Magnetic field frequency range [Hz]:	1-100
Magnetic field intensity [mT]:	0-20
Pulse/break duration [s]:	0,5-8
Power supply [V/Hz/W]:	230/50/400
Weight of control unit [kg]:	6,2
Dimensions of control unit (H x W x L) [cm]:	14,2 x 36,4 x 33,5
Shapes of magnetic field modulation:	sine, rectangle, triangle – unipolar and bipolar
Magnetic field applicators (optional):	AS-200K, AS-315K, AS-600K, AS-200N, AS-315N, AS-600N, AP-100, APP-100, AP-ID, APP-2D
Mobile trolleys for applicators (optional):	S-200N, S-315N
L-6 couch for applicators (optional):	AS-600K, AS-600N



MULTITRONIC MT-4-MED INMED



★ FEATURES

- Modern device for two-channel electrotherapy and laser therapy.
- Two simultaneous treatments option.
- Colour touchscreen display (4.3").
- User-friendly touchscreen and button controls.
- Ergonomic and lightweight probes.
- Ready to use preset treatment programs for popular treatments.
- Easy to store user-defined programs thanks to on-screen keyboard.
- Independent adjustment of treatment settings.
- Fan control feature to minimize noise and power consumption.
- Treatment duration and number of performed treatments counter.
- Can be used as portable device.

TECHNICAL DATA

Power supply:	~230V 10%, 50Hz, 70VA
Electrical protection class:	I, type BF
Dimensions [mm]:	335 x 270 x 125 mm
Weight (of control unit):	3.1 kg
Laser class:	3B
Laser treatment timer:	1 s -99 min



DEEP MUSCLES STIMULATION ZEPU



ZEPU-K5000A

Deep muscle stimulator

★ FEATURES

- Vibration and rebounding technology.
- Adjustable frequency technology.
- Several hand pieces for option.
- Small, light weight and portable with handbag.

ZEPU-K5000A SPECIFICATIONS

Size (mm): 266×50×146(L×W×H)

Weight: 2.5kgs



ZEPU-L100

Cryo therapy apparatus

★ FEATURES

- Fast cooling and significant effect.
- Pinpoint the accurate treatment area.
- Special used for acute sports injury.

ZEPU-L100 SPECIFICATIONS

Size (mm): 320×585×870(L×W×H)

Weight: 55kgs



DYNAMIC INTERFERENCE ELECTRIC ZEPU



ZEPU-G60

Stereoscopic dynamic interference electric therapeutic apparatus

★ FEATURES

- Three-way output, 6-electrode three-dimensional space stimulation effect.
- Independent dual channel.
- Adsorption electrode.

ZEPU-G60 SPECIFICATIONS

Size (mm): 950x338x458(L×W×H)

Weight: 25kgs



ZEPU-K2

Upper and lower limb air pressure therapy apparatus

★ FEATURES

- 6-chamber airbags, stacked design, effectively prevent back-flow of body fluids.
- 10 working modes, 1 standard mode and 10 extended modes.
- The host can connect 2 inflatable airbags at the same time.
- The host adopts special shock absorbing device with low noise.

ZEPU-K2 SPECIFICATIONS

Size (mm): 360×270×175(L×W×H)

Weight: 12kgs



INFRARED ZEPU

ZEPU-GL575/750

Infrared Therapeutic Apparatus

★ FEATURES

- Mobile floor cabinet and lifting treatment head design.
- Electric lifting easy to set the height.
- Color LCD touch screen.
- User-defined treatment plan and store it.



ZEPU-GL575 SPECIFICATIONS

Size (mm): 560 x 520 x (1200-1600) (L×W×H)

Weight: 30kgs



ZEPU-GL750

ZEPU-GL750 SPECIFICATIONS

Size (mm): 500×500× (1200-1600) (L×W×H)

Weight: 65kgs



MEDICAL LASER THERAPY

ELECTRONIC ZEPU



★ PRODUCT PRINCIPLE

When semiconductor laser acts on the body, it will produce a biological stimulating effect, which can promote cell regeneration, accelerate tissue repair, improve blood circulation, reduce inflammation, relieve pain, reduce edema, and regulate the body's immune function. It is not only suitable for superficial skin lesions, but also for deep tissue lesions, showing a therapeutic effect from shallow to deep.

★ FEATURES

- Combining treatment at selected spots with that in the entire area
Double treatment head, square treatment head has large treatment area, pen type treatment head has focusing function and strong penetrability.
- Dual-channel output
Two independent outputs, namely, A and B, which can be independently controlled
- Streamlined operations
Large LCD touch screen, humanized operation interface; cylinder hydraulic support can be accurately adjusted in all directions, simple and flexible.
- Stable performance
GaAlAs laser and professional circuit design have superior performance and high stability.

SPECIFICATIONS

ZYS-808A	ZYS-808	ZYS-808B
Single output	Dual-channel output	Dual-channel output
Square treatment head	Square therapeutic head+pen therapeutic head	2 square treatment heads



SHOCK WAVE THERAPY SYSTEM ZEPU



ZEPU-C2

Extracorporeal shock wave therapy system

★ FEATURES

- Independent dual channel.
- Double channels and double handles, 6 treatment handpieces and 2 massage handpieces for options.
- Provide automatic mode and manual selection mode.

ZEPU-C2 SPECIFICATIONS

Size (mm): 400x320x160(LxWxH)

Weight: 8kgs



ZEPU-D50

Low frequency electronic pulse therapy apparatus

★ FEATURES

- Promote local blood circulation in the human body, delay muscle atrophy.
- Using vacuum adsorption pressure, the adsorption pressure can be adjusted.

ZEPU-D50 SPECIFICATIONS

Size (mm): 410x345x180(LxWxH)

Weight: 11.2kgs



ZEPU-P20

High-frequency chest wall oscillation sputum meter

★ FEATURES

- Color high-definition LCD display, navigation operation guide.
- Provide 3 sizes of inflatable vest to choose from.
- Pneumatic emergency stop safety switch.
- Optional automatic detection and control modules.
- Configure the smart bluetooth wireless connection port.

ZEPU-P20 SPECIFICATIONS

Size (mm): 350x380x240 (LXWXH)

Weight: 15kgs

ZEPU-K5000A

Deep muscle stimulator



ZEPU-K5000A

Deep muscle stimulator

★ FEATURES

- Vibration and rebounding technology.
- Adjustable frequency technology.
- Several hand pieces for option.
- Small, light weight and portable with handbag.

ZEPU-K5000A SPECIFICATIONS

Size (mm): 266x50x146(LxWxH)

Weight: 2.5kgs



ZEPU-L100

Cryo therapy apparatus

★ FEATURES

- Fast cooling and significant effect.
- Pinpoint the accurate treatment area.
- Special used for acute sports injury.

ZEPU-L100 SPECIFICATIONS

Size (mm): 320x585x870(LxWxH)

Weight: 55kgs



ZEPU-G60

Stereoscopic dynamic interference electric therapeutic apparatus



ZEPU-G60

Stereoscopic dynamic interference electric therapeutic apparatus

★ FEATURES

- Three-way output, 6-electrode three-dimensional space stimulation effect.
- Independent dual channel.
- Adsorption electrode.

ZEPU-G60 SPECIFICATIONS

Size (mm): 950x338x458(LxWxH)

Weight: 25kgs



ZEPU-K2

Upper and lower limb air pressure therapy apparatus

★ FEATURES

- 6-chamber airbags, stacked design, effectively prevent back-flow of body fluids.
- 10 working modes, 1 standard mode and 10 extended modes. The host can connect 2 inflatable airbags at the same time. The host adopts special shock absorbing device with low noise.

ZEPU-K2 SPECIFICATIONS

Size (mm): 360x270x175(LxWxH)

Weight: 12kgs



ZEPU-C2,ZEPU-D50 AND ZEPU-P20



ZEPU-C2

Extracorporeal shock wave therapy system

★ FEATURES

- Independent dual channel.
- Double channels and double handles, 6 treatment handpieces and 2 massage handpieces for options.
- Provide automatic mode and manual selection mode.

ZEPU-ZD1 SPECIFICATIONS

Size (mm): 400x320x160 (LxWxH)

Weight: 8kgs



ZEPU-D50

Low frequency electronic pulse therapy apparatus

★ FEATURES

- Promote local blood circulation in the human body, delay muscle atrophy.
- Using vacuum adsorption pressure, the adsorption pressure can be adjusted.

ZEPU-D50 SPECIFICATIONS

Size (mm): 410x345x180(LxWxH)

Weight: 11.2kgs



ZEPU-P20

High-frequency chest wall oscillation sputum meter

★ FEATURES

- Color high-definition LCD display, navigation operation guide.
- Provide 3 sizes of inflatable vest to choose from. Pneumatic emergency stop safety switch.
- Optional automatic detection and control modules. Configure the smart bluetooth wireless connection port.

ZEPU-P20 SPECIFICATIONS

Size (mm): 350x380x240(LXWXH)

Weight: 15kgs

ZEPU-ZD1,2,3,5,6

VERTICAL WAVE VIBRATION SERIES

★ FEATURES

- Vertical sound wave vibration.
- No working noisy.
- Accelerated speed could be controlled.
- Adjustable vibration intensity.
- More accurate and high reliability.



ZEPU-ZD1

Vertical wave vibration

ZEPU-ZD1 SPECIFICATIONS

Size (mm):	840x860x1480 (LxWxH)
Weight:	87kgs



ZEPU-ZD2

Vertical wave vibration

ZEPU-ZD2 SPECIFICATIONS

Size (mm):	840x860x1480 (LxWxH)
Weight:	87kgs



ZEPU-ZD3

ZEPU-ZD3 SPECIFICATIONS

Size (mm):	355x1530x1070 (LxWxH)
Weight:	240kgs



ZEPU-ZD5

ZEPU-ZD5 SPECIFICATIONS

Size (mm):	1980x800x650 (LxWxH)
Weight:	145kgs



ZEPU-ZD6

ZEPU-ZD5 SPECIFICATIONS

Size (mm):	870x770x1070 (LxWxH)
Weight:	72kgs

WHY CHOOSE IPEX HEALTHCARE?



Complete
Rehabilitation
Solutions



Internationally
Recognized
Manufacturers



Clinical Training
& Technical
Support



Preventive
Maintenance &
After-sales Service



Equipment Supply,
Installation &
Commissioning



Nationwide
Sales & Technical
Coverage



Postal Address:

IPEX Healthcare Limited,
P. O. Box 42502, Dar Es Salaam, Tanzania.

Physical Address:

Plot No. 08, Bahati Street, Temeke Area,
Dar es Salaam, Tanzania

Tel: +255 712 897 082 / +255 745 226 582

Email: info@ipex.co.tz **Web:** www.ipex.co.tz