



Rehabilitation Assessment & Diagnostics

PRODUCT CATALOGUE



Scan to Download
the Catalogue

YOUR RELIABLE & TRUSTED PARTNER

**In Physiotherapy, Rehabilitation
& Occupational Health Solutions**



ABOUT US

IPEX Healthcare Limited is a Tanzanian based and Tanzanian registered company dealing with Importation, and Supply of Hospital consumables, Medical, Surgical and Laboratory equipment, supplies, devices and other related products and services in the medical and healthcare sector.

STORAGE CART FERROX

★ FEATURES

(1) Art. 260610 – Storage CART **

- natural wood
- for clubs, hoops and balls (etc.)
- dimensions 150 x 60 x 104 h cm

Version shown with optional accessories Art. 260615

(2) Art. 260615 – CART EQUIPMENT KIT ** equipped with:


- clubs, 36 cm
- clubs, 45 cm
- sand bags, 38 x 20 cm
- solid wood "Vayer" bricks
- plastic blocks for creating hurdle courses
- coloured nylon ropes
- ropes with handles and thickening
- balls for rhythmic gymnastics 130 g, 16 cm
- balls for rhythmic gymnastics 280 g, 18 cm
- balls for rhythmic gymnastics 420 g, 19 cm
- PVC gymnastic ring, 60 cm
- PVC gymnastic ring, 70 cm
- PVC gymnastic ring, 80 cm
- Bobath ball, 55 cm
- Bobath ball, 75



PSYCHOMOTOR TRAINING FISIOTECH

HOOPS

Colourful nylon hoops, flat section. Used for psychomotor courses along with cones with holes (code 1416031) and sticks (code 1416061).




SPECIFICATIONS

MODELS	COLOURS	CODE
Ø 60 cm	assorted colours	1416011
Ø 70 cm	assorted colours	1416012
Ø 80 cm	assorted colours	1416013

PVC STICKS SET – CODE 1416061

Set of four colourful gymnastics PVC sticks. Used for psychomotor courses along with cones with holes (code 1416031) and hoops (codes 1416011–12–13).



SPECIFICATIONS

MODELS	COLOURS	CODE
Ø 3 cm	assorted colours	1416061

CONE WITH HOLES – CODE 1416031

Cone with holes available in four random colours (red, yellow, green, blue). Split vertex to insert PVC sticks (code 1416061) and hoops (codes 1416011–12–13) and three symmetrical lateral holes to insert PVC sticks in order to create courses.



SPECIFICATIONS

MODELS	COLOURS	CODE
30 cm H	assorted colours	1416011

BALANCE ASSESSMENT AND REHABILITATION-NEUROSOFT

★ ASSESSMENT

- ONE SENSOR FOR COMPREHENSIVE BALANCE ANALYSIS
- COM ANALYSIS BASED ON TIME AND FREQUENCY
- SCREENING AND ADVANCED EXAMINATIONS
- ALL STANDARD CLINICAL PROTOCOLS (LOS TEST, ROMBERG'S TEST, ETC.)



GAIT ANALYSIS AND TRAINING WITH BIOFEEDBACK -NEUROSOFT

ASSESSMENT


- GAIT ANALYSIS FROM SIMPLE SCREENING TO A COMPREHENSIVE EXAM
- UP TO 25 GAIT PARAMETERS (TEMPORAL, SPATIAL, 3D KINEMATIC) AND DYNAMIC SURFACE EMG FROM UP TO 24 MUSCLES
- INTERFACE OPTIMIZED FOR CLINICAL USE



10 min turnaround between a patient to come and report generation



Signal acquisition is not affected by the surrounding metal constructions




JOINT MOVEMENT ANALYSIS-NEUROSOFT

ASSESSMENT

- ANALYSIS USING IMU SENSORS
- SYNCHRONIZED EMG AND JOINT MOVEMENT ANALYSIS
- DATA EXCHANGE VIA WI-FI



NEUROSENS – NEUROSOF



NEUROSENS



All-in-one. Acquisition of gait parameters and EMG is done with one IMU sensor



The data exchange with computer software is done via Wi-Fi

REHABILITATION

- TRAINING OF COORDINATION AND MOVEMENTS IN JOINT AND SPINE
- EARLY REHABILITATION POSSIBILITIES
- MANUAL OR AUTOMATIC TRAINING MODES

ZEPU-SG2-3

Digital OT assessment and training system

ZEPU-SG2/3

Digital OT assessment and training system

★ FEATURES

SG2: Provide more than 50 interactive training modes
Able to perform finger fine motor control training, upper limb motor function training and evaluation recording.
SG3: A variety of accessories are available to evaluate and train fine movements of upper limbs and hands.



ZEPU-SG2



ZEPU-SG3

ZEPU-RX1

Multi-joint flexible comprehensive rehabilitation trainer

★ FEATURES

Multiple movements of the wrist joint can be trained, palmar flexion, dorsiflexion, radial deviation, ulnar deviation, pronation, and supination.



ZEPU-SG2/3 SPECIFICATIONS

Size (mm):

1070x800x(1050-1350) (LxWxH)

Weight:

61kgs/75kgs

WHY CHOOSE IPEX HEALTHCARE?



Complete
Rehabilitation
Solutions



Internationally
Recognized
Manufacturers



Clinical Training
& Technical
Support



Preventive
Maintenance &
After-sales Service



Equipment Supply,
Installation &
Commissioning



Nationwide
Sales & Technical
Coverage



Postal Address:

IPEX Healthcare Limited,
P. O. Box 42502, Dar Es Salaam, Tanzania.

Physical Address:

Plot No. 08, Bahati Street, Temeke Area,
Dar es Salaam, Tanzania

Tel: +255 712 897 082 / +255 745 226 582

Email: info@ipex.co.tz **Web:** www.ipex.co.tz